

国内発売記念キャンペーン 15% off

大学・大学院・研究機関 様

☞ 平成 28 年 2 月 8 日より Magritek 社 卓上型 NMR 冷媒不要 販売開始！
平成 29 年 3 月 31 日まで実施します。

Spinsolve は
卓上型 NMR ですから
場所を問いません。実
験室・ラボ・ドラフト
チャンバー等に設置
できます。

学生実験・合成リアル
タイムモニター・スク
リーニングなどに最
適な NMR となりま
す。ヘリウムガス等は
不要ですのでランニ
ングコストも安価で
予算にも優しい NMR
です。

機種：下記 4 種です。

1H 専用

1H, 19F 2D NMR

1H, 19F, 13C 2D NMR

1H, 19F, 31P 2D NMR

主な仕様：

周波数 42.5MHz(1H)

寸法：L58xW43xH40cm

重量：約 55Kg

磁石：永久磁石

High-resolution NMR spectroscopy to monitor reactions under the fume hood

Spinsolve™
The fastest compact NMR spectrometer

magritek

www.magritek.com/spinsolve

FLOW SETUP UNDER THE FUME HOOD

In this work we tested the performance of the benchtop system under the fume hood of a chemistry lab. It was used to monitor the conversion of two different reactions in real time. Taking advantage of the clear bore of the magnet a flow cell was mounted in the magnet and a closed flow setup was used to pump the reacting mixtures from the reactor through the magnet and back. A peristaltic pump sets the flow rate with high accuracy. Before running the reactions the effect of flow was investigated to determine the maximum flow rate that does not introduce any line broadening in the spectra.



Figure 1. Experimental set up. (a) NMR magnet; (b) NMR flow cell; (c) Peristaltic pump to provide the continuous flow of the reaction mixture through the magnet bore. The NMR spectrometer is outside the fume hood and cannot be seen in the pictures.

MIXING MONITORING AND CONCENTRATIONS DETERMINATION

The benchtop high-resolution NMR setup was tested first to measure the concentration ratio in mixtures with different dilution degrees. The solvents were pumped at variable flow-rates and the NMR signal was recorded on-line every 10 seconds while the sample passed through the magnet (Fig. 2a) through a flow cell 4 mm in diameter and 20 mm long. Figure 2b shows that after varying the flow rate in the pumps a short transitory of about one minute was enough to reach the new signal value proportional to the known concentration. Figure 2c plots the correlation between flow rate (pump concentration) vs. the NMR signal. The highly linear correlation measure for both solvents demonstrate the accuracy of the system to quantify concentration even under continuous flow.




Figure 2. (a) Sketch of the flow-NMR setup using two pumps to mix different solutions. (b) Concentrations of solubates A and B determined from the NMR spectra as a function of time. (c) Correlation of the concentration of solutions A and B obtained from the NMR spectra and the ones calculated by the flow rates.

References

- J. Peña, F. Casanova, and B. Blumrich, *Science*, 2007, 315, 1110-1112.
- E. Daniele, J. Peña, B. Blumrich, and F. Casanova, *Angew. Chem. Int. Ed.*, 2010, 49, 4132-4135.
- S. K. K. sey, E. Daniele, B. Blumrich, and F. Casanova, *Phys. Chem. Chem. Phys.*, 2011, 13, 13172.

REAL TIME MONITORING OF BIODIESEL PRODUCTION

To obtain a high quality fuel and to achieve the requirements for commercialization, it is important to monitor the transesterification process to ensure low levels of mono- and diacylglycerols resulting from incomplete reactions. In this work we show the potential of the present system to monitor the progress of the catalyzed transesterification reaction of vegetable oils that produces biodiesel. To perform a previous test a transesterification reaction of rape oil was carried out in a 500 mL three-necked glass flask and continuously analyzed with the benchtop high-resolution NMR system by a continuous flow of the reaction mixture through the magnet bore. Figure 3 shows the experimental set up.




Figure 3. (a) Transesterification reaction equation. In our experiment 200 mL of rape seeds oil, 75 mL of methanol (0:1 M ratio) and 1.35 g of NaOH (0.5% m/m related to oil) were used. (b) Top: Spectra of pure rape oil and pure biodiesel. Bottom: Spectra monitored during 60 minutes while the reaction proceeded. (c) Zoomed spectra of the region of interest.

CATALYTIC REDUCTION OF ACETOPHENONE BY ISOPROPANOL

Figure 4 shows another reaction that can be monitored by the benchtop system. The reduction of Acetophenone was carried out on a mixture with very small concentration (0.5 Molar). This can be appreciated on the full scale spectrum, where mainly the isopropanol lines can be identified. However, zooming in the spectra acquired as a function of time, the concentration of products and reactants can be clearly quantified. Besides demonstrating that the NMR system offers a remarkable performance to study a process of high industrial relevance, these results prove that very small concentrations can be quantified with high accuracy even under continuous flow conditions.




Figure 4. (a) Sketch of the reduction induced by the catalyst. (b) Full scale and zoom of the spectra at the beginning and at the end of the monitored time windows. (c) Time dependence of the lines in the aromatic region. Each spectrum is an average of 8 scans. The lines can be integrated to quantify each component (b).

☎ **Spinsolve** についてのお問合せ・ご質問は下記へお問合せ願います。



info@asahilab.co.jp



弊社は下記の製品を輸入販売しておりますので
合わせましてご案内申し上げます。
何卒宜しくお願い申し上げます。

Remote control ラボジャッキ



卓上型 NMR Spinsolve

販売代理店募集中！

- ☑ **Bochem**(ドイツ) ラボジャッキ・ボスヘッド・ブンセンバーナー・ステンレス家具他
- ☑ **FutureChemistry** (オランダ) マイクロフローリアクター合成装置
- ☑ **Spinsolve** (ニュージーランド・ドイツ) **卓上型 NMR 冷媒不要!**

詳細は下記へご連絡下さい。

info@asahilab.co.jp

078-335-8613 Tel

078-335-8614 Fax

<http://www.asahilab.co.jp/>

〒650-0044 神戸市中央区東川崎町 1-8-4-801



マイクロフロー合成装置
Integrated with IR

For future science

輸入販売元：



株式会社朝日ラボ交易